

Determining Fatty Acid – Membrane Partition Coefficients

Synopsis

ADIFAB can be used to determine the partition coefficient of a fatty acid between a membrane phase and aqueous solution. Simply add fatty acid to a cuvette containing ADIFAB and a membrane and measure the fluorescence ratio (505/432 upon excitation at 386 nm).

Procedure

For details on measuring the ADIFAB ratio and calculating [FFA] see [Determining the ADIFAB Ratio](#). To determine R_0 , add ADIFAB and a membrane of known lipid concentration to a cuvette containing measuring buffer (20 mM HEPES, 140 nM NaCl, 5 mM KCl, 1 mM Na_2HPO_4 , at pH 7.4), and measure the fluorescence ratio (505/432 nm). Titrate the solution with fatty acid aliquots of known concentration and measure the R value after each addition—be sure to wait at least 5-10 minutes for equilibrium before measuring R. For each R measured, calculate the partition coefficient, K_p , using Eq. (1):

$$K_p = \frac{\frac{V_a}{V_m} \cdot ([\text{FA}]_{\text{total}} - [\text{FFA}])}{[\text{FFA}]} \quad (1)$$

where $[\text{FA}]_{\text{total}}$ is the *total* fatty acid concentration in the cuvette after each addition, $[\text{FFA}]$ is the *free* fatty acid concentration, and V_a and V_m are the volumes of the aqueous and membrane phases, respectively. (Note: to accurately report $[\text{FA}]_{\text{total}}$, measure the concentration of the fatty acid stock according to [Determining the Concentration of Fatty Acid in an Aqueous Solution](#).) With sufficient membrane present, no correction for wall binding to the cuvette walls is necessary because for typical conditions >95% of the fatty acid will be bound to the membrane; very little will be free, bound to ADIFAB or bound to the walls.

Example

In a cuvette containing 2 ml buffer (20 mM HEPES, 140 nM NaCl, 5 mM KCl, 1 mM Na_2HPO_4 , at pH 7.4 and 37°C), 100 μM egg phosphatidylcholine vesicles (EPC) and 0.2 μM ADIFAB, the R_0 was measured and found to be 0.279. 3 μM of sodium palmitate was added, and after waiting 10 minutes for equilibrium, the R value was measured and found to be 0.335. Using the R value to determine that $[\text{FFA}] = 88.6$ nM and substituting $V_a/V_m = 10000$ and $[\text{FA}]_{\text{total}} = 3$ μM into Eq. (1), K_p was calculated to be 3.29×10^5 . Additional aliquots of sodium palmitate were added and after each addition R values were measured and K_p values were calculated. The average value of K_p for the titration was 3.47×10^5 . The complete set of R, $[\text{FFA}]$ and K_p values are listed in Table 1.

Table 1. Titration data from palmitate and EPC vesicles K_p determination.

[FA] _{total} (uM)	Measured R Value	[FFA] (nM)	$K_p \times 10^{-5}$
0	0.279	0	
3	0.427	89	3.29
5	0.516	143	3.40
7	0.606	199	3.42
10	0.718	270	3.60
15	0.920	402	3.63
20	1.150	558	3.48
25	1.345	696	3.49
		average =	3.47